


**Tubos para recoger muestras de gas, tubo DURAN<sup>®</sup>****Lenz****Tubos para recoger muestras de gas, tubo DURAN<sup>®</sup>**

Fabricados a partir de tubo DURAN<sup>®</sup>. Disponibles con o sin boquilla (con rosca GL y septum), con dos llaves de paso unidireccionales NS14 (con dispositivo de retención), conexiones: tubo capilar (diámetro exterior 8mm, diámetro interior 2,7mm), oliva de vidrio (diámetro exterior 8mm), con etiqueta de inscripción. También disponibles con revestimiento PUR.



Code	Description	Packaging
LLG06243069	Gas sampling tube 150mlwith extraction	1 pz.
LLG07623930	Gas sampling tube 1000mlwith extraction	1 pz.
LLG09110215	Gas Collecting Tubes without Nozzle, Ca	1 pz.
LLG09110220	Gas Collecting Tubes without Nozzle, Ca	1 pz.
LLG09110225	Gas Collecting Tubes without Nozzle, Ca	1 pz.
LLG09110230	Gas Collecting Tubes without Nozzle, Ca	1 pz.
LLG09110236	Gas collecting tubes,DURAN®,with stopco	1 pz.
LLG09110237	Gas collecting tubes,DURAN®,with stopco	1 pz.



**CARLO ERBA Reagents SA**  
 Calle Filadors 35,  
 6º Planta 5º Puerta  
 08208 Sabadell (BCN)  
 Tel. +34 93 693 37 35  
[www.carloerbareagents.com](http://www.carloerbareagents.com)



**CARLO ERBA**  
*Reagents operates with  
 a Certified Quality  
 Management System*

