



Tapones roscados y accesorios para bidones, HDPE

Kautex

Tapones roscados y accesorios para bidones, HDPE



Code	Description	Packaging
LLG09140040	Screw cap, HDPE with detachable ring, r	1 pz.
LLG09140042	Screw cap, HDPE with detachable ring, r	1 pz.
LLG09140046	Screw cap black, PE-HD, with tearing ri	1 pz.
LLG09140048	Screw cap, black, with tear-off ringfor	1 pz.
LLG09140050	Screw cap, black, with tear-off ringfor	1 pz.
LLG09140051	Screw cap, black, with tear-off ringfor	1 pz.
LLG09140054	Closure with dispensing pump, PE-HDfor	1 pz.
LLG09140061	Closure for 5 and 10 ltr type	1 pz.
LLG09140062	Screw cap with drain cockHDPE for 2, 2,	1 pz.
LLG09140063	Closure for 20 and 30 ltr type	1 pz.
LLG09140065	Screw cap PE-HD with dosing pumpfor 10l	1 pz.
LLG09140068	Screw cap PE-HD, blackwith tear-off rin	1 pz.
LLG09140070	Screw-cap, PE-HD, with dosing pumpfor 2	1 pz.
LLG09140071	Screw cap HDPE, with dosing pumpfor 25	1 pz.
LLG09140072	Srew cap HDPE with dosing pumpfor 30L	1 pz.
LLG09140073	Key, HDPEDIN45 black,for 2, 2,5, 3 and	1 pz.
LLG09140108	Jerrycan key, HDPEfor 5,10,20,25 and 30	1 pz.





CARLO ERBA Reagents SA
Calle Filadors 35,
6º Planta 5º Puerta
08208 Sabadell (BCN)
Tel. +34 93 693 37 35
www.carloerbareagents.com



CARLO ERBA
*Reagents operates with
a Certified Quality
Management System*

