

**Baños termostáticos con refrigeración a -55°C****Huber****Baños termostáticos con refrigeración a -55°C**

Forma compacta y elevada potencia de refrigeración a temperaturas bajas.

Opcional: refrigerante natural disponible bajo pedido

Potencia de calefacción

CC-505/CC-505wl:

CC-508 hasta CC-525w:

1,5 kW

3,0 kW

Datos de la bomba CC-505 a CC-525w

Presión max.:

Succión max.:

25 l/min. / 0,7 bar

18,5 l/min. / 0,4 bar

Requisitos de alimentación

CC-505 a CC-508w:

CC-510 a CC-525w:

230 V / 50 Hz

400 V / 50 Hz

Estabilidad de la temperatura según DIN 12876: +/-0,02K



Code	Description	Packaging
LLG09857401	Refrigerated heating bath circulator CC	1 pz.
LLG09859241	Refrigerated heating bath circulator CC	1 pz.
LLG09859242	Refrigerated heating bath circulator CC	1 pz.
LLG09859245	Refrigerated heating bath circulator CC	1 pz.
LLG09859246	Refrigerated heating bath circulator CC	1 pz.
LLG09859247	Refrigerated heating bath circulator CC	1 pz.
LLG09859248	Refrigerated heating bath circulator CC	1 pz.
LLG09859249	Refrigerated heating bath circulator CC	1 pz.
LLG09859256	Cold-Warmth-Bath-Circulation thermostat	1 pz.



CARLO ERBA Reagents SA
 Calle Filadors 35,
 6º Planta 5º Puerta
 08208 Sabadell (BCN)
 Tel. +34 93 693 37 35
www.carloerbareagents.com



CARLO ERBA
*Reagents operates with
 a Certified Quality
 Management System*

