

**Baños termostáticos con refrigeración a -55°C****Huber****Baños termostáticos con refrigeración a -55°C**

Forma compacta y elevada potencia de refrigeración a temperaturas bajas.

Opcional: refrigerante natural disponible bajo pedido

Potencia de calefacción

CC-505/CC-505wl:

CC-508 hasta CC-525w:

1,5 kW

3,0 kW

Datos de la bomba CC-505 a CC-525w

Presión max.:

Succión max.:

25 l/min. / 0,7 bar

18,5 l/min. / 0,4 bar

Requisitos de alimentación

CC-505 a CC-508w:

CC-510 a CC-525w:

230 V / 50 Hz

400 V / 50 Hz

Estabilidad de la temperatura según DIN 12876: +/-0,02K



| Code | Description | Packaging |
|-------------|-----------------------------------------|-----------|
| LLG09857401 | Refrigerated heating bath circulator CC | 1 pz. |
| LLG09859241 | Refrigerated heating bath circulator CC | 1 pz. |
| LLG09859242 | Refrigerated heating bath circulator CC | 1 pz. |
| LLG09859245 | Refrigerated heating bath circulator CC | 1 pz. |
| LLG09859246 | Refrigerated heating bath circulator CC | 1 pz. |
| LLG09859247 | Refrigerated heating bath circulator CC | 1 pz. |
| LLG09859248 | Refrigerated heating bath circulator CC | 1 pz. |
| LLG09859249 | Refrigerated heating bath circulator CC | 1 pz. |
| LLG09859256 | Cold-Warmth-Bath-Circulation thermostat | 1 pz. |



CARLO ERBA Reagents SA
 Calle Filadors 35,
 6º Planta 5ª Puerta
 08208 Sabadell (BCN)
 Tel. +34 93 693 37 35
www.carloerbareagents.com



CARLO ERBA
*Reagents operates with
 a Certified Quality
 Management System*

