

**Termostatos de frío de control compatibles hasta -90 °C****Huber****Termostatos de frío de control compatibles hasta -90 °C**

El modelo CC-805 es una alternativa económica para bajas temperaturas cuando se necesita poca potencia.

Opción: refrigerante natural por encargo

Datos de la bomba CC-805 hasta CC-906w:

Presión máx.:

Succ. máx.: 25 l/min / 0,7 bar

18,5 l/min / 0,4 bar

Conexión de red:

CC-805 / CC-902: 400 V / 50 Hz

230 V / 50 Hz

Estabilidad de temperatura según DIN 12876: 0,02 K



Code	Description	Packaging
LLG09857400	Refrigerated heating bath circulator CC	1 pz.
LLG09859250	Refrigerated heating bath circulator CC	1 pz.
LLG09859251	Refrigerated heating bath circulator CC	1 pz.
LLG09859252	Refrigerated heating bath circulator CC	1 pz.
LLG09859253	Refrigerated heating bath circulator CC	1 pz.
LLG09859254	Refrigerated heating bath circulator CC	1 pz.
LLG09859255	Refrigerated heating bath circulator CC	1 pz.



CARLO ERBA Reagents SA
Calle Filadors 35,
6º Planta 5ª Puerta
08208 Sabadell (BCN)
Tel. +34 93 693 37 35
www.carloerbareagents.com



CARLO ERBA
*Reagents operates with
a Certified Quality
Management System*

