

**Termostatos de frío de control compatibles hasta -90 °C****Huber****Termostatos de frío de control compatibles hasta -90 °C**

El modelo CC-805 es una alternativa económica para bajas temperaturas cuando se necesita poca potencia.

Opción: refrigerante natural por encargo

Datos de la bomba CC-805 hasta CC-906w:

Presión máx.:

Succ. máx.: 25 l/min / 0,7 bar

18,5 l/min / 0,4 bar

Conexión de red:

CC-805 / CC-902: 400 V / 50 Hz

230 V / 50 Hz

Estabilidad de temperatura según DIN 12876: 0,02 K



| Code | Description | Packaging |
|-------------|---|-----------|
| LLG09857400 | Refrigerated heating bath circulator CC | 1 pz. |
| LLG09859250 | Refrigerated heating bath circulator CC | 1 pz. |
| LLG09859251 | Refrigerated heating bath circulator CC | 1 pz. |
| LLG09859252 | Refrigerated heating bath circulator CC | 1 pz. |
| LLG09859253 | Refrigerated heating bath circulator CC | 1 pz. |
| LLG09859254 | Refrigerated heating bath circulator CC | 1 pz. |
| LLG09859255 | Refrigerated heating bath circulator CC | 1 pz. |



CARLO ERBA Reagents SA
 Calle Filadors 35,
 6º Planta 5ª Puerta
 08208 Sabadell (BCN)
 Tel. +34 93 693 37 35
www.carloerbareagents.com



CARLO ERBA
*Reagents operates with
 a Certified Quality
 Management System*

