



**Areómetro, para lechada de cal  $\text{Ca(OH)}_2$**

**Areómetro, para lechada de cal  $\text{Ca(OH)}_2$**

Para lechada de cal  $\text{Ca(OH)}_2$ .

Sin termómetro.

Temperatura de referencia +20 °C.



Code	Description	Packaging
LLG09004640	Hydrometer,for Ca(OH) <sup>2</sup> ,without thermome	1 pz.



**CARLO ERBA Reagents SA**  
Calle Filadors 35,  
6º Planta 5º Puerta  
08208 Sabadell (BCN)  
Tel. +34 93 693 37 35  
[www.carloerbareagents.com](http://www.carloerbareagents.com)



**CARLO ERBA**  
*Reagents operates with  
a Certified Quality  
Management System*

