



Fascos graduados de cuello estrecho

Fascos graduados de cuello estrecho

Material: PE. Ideales para agua y aceites: no tienen bordes que impidan verter los líquidos.

Graduados y, por tanto, útiles para mediciones aproximadas. Hombros perforados para precintarse el frasco y evitar la contaminación de las muestras durante el transporte. Según DIN 13316 y 168. Adecuado para alimentos.

Pida los tapones roscados por separado.



Code	Description	Packaging
LLG09102433	Narrow neck bottle 50ml, PEG L 18, gradu	1 pz.
LLG09102434	Narrow neck bottle 100ml, PEG L 18, grad	1 pz.
LLG09102435	Narrow neck bottle 250ml, PEG L 25, grad	1 pz.
LLG09102436	Narrow neck bottle 500ml, PEG L 25, grad	1 pz.
LLG09102437	Narrow neck bottle 1000ml, PEG L 32, gra	1 pz.
LLG09102438	Narrow neck bottle 2000ml, PEG L 32, gra	1 pz.



CARLO ERBA Reagents SA
 Calle Filadors 35,
 6º Planta 5ª Puerta
 08208 Sabadell (BCN)
 Tel. +34 93 693 37 35
www.carloerbareagents.com



CARLO ERBA
*Reagents operates with
 a Certified Quality
 Management System*

