



**Tapones roscados, HDPE**

**Kartell**

**Tapones roscados, HDPE**

Tapón roscado estándar para frascos DIN168.



Code	Description	Packaging
LLG09102439	Screw cap GL 18, PEfor 9102433 and 9102	1 pz.
LLG09102440	Screw cap GL 25, PEfor 9102435 and 9102	1 pz.



**CARLO ERBA Reagents SA**  
 Calle Filadors 35,  
 6º Planta 5º Puerta  
 08208 Sabadell (BCN)  
 Tel. +34 93 693 37 35  
[www.carloerbareagents.com](http://www.carloerbareagents.com)



**CARLO ERBA**  
*Reagents operates with  
 a Certified Quality  
 Management System*

