



Tapones roscados, HDPE

Kartell

Tapones roscados, HDPE

Tapón roscado estándar para frascos DIN168.



Code	Description	Packaging
LLG09102439	Screw cap GL 18, PEfor 9102433 and 9102	1 pz.
LLG09102440	Screw cap GL 25, PEfor 9102435 and 9102	1 pz.



CARLO ERBA Reagents SA
Calle Filadors 35,
6º Planta 5º Puerta
08208 Sabadell (BCN)
Tel. +34 93 693 37 35
www.carloerbareagents.com



CARLO ERBA
*Reagents operates with
a Certified Quality
Management System*

