

**Vasos de precipitados graduados, forma baja****Kartell****Vasos de precipitados graduados, forma baja**

Polipropileno. Según ISO 7056 - 1981 (E). Polipropileno de nueva generación que proporciona excelente translucidez. Autoclavable a 121°C. Graduación en relieve hecha en el mismo molde para garantizar la reproducibilidad. Excelente resistencia química.



Code	Description	Packaging
LLG09013556	Beaker 50ml, PPgraduated, low formautoc	1 pz.
LLG09013557	Beaker 100ml, PPgraduated, low formauto	1 pz.
LLG09013558	Beaker 250ml, PPgraduated, low formauto	1 pz.
LLG09013559	Beaker 500ml, PPgraduated, low formauto	1 pz.
LLG09013560	Beaker, 1000ml, PPgraduated, low formau	1 pz.
LLG09013561	Beaker 2000ml, PPgraduated, low formaut	1 pz.
LLG09013563	Beaker 5000ml, PPgraduated, low formaut	1 pz.



**CARLO ERBA Reagents SA**  
 Calle Filadors 35,  
 6º Planta 5º Puerta  
 08208 Sabadell (BCN)  
 Tel. +34 93 693 37 35  
[www.carloerbareagents.com](http://www.carloerbareagents.com)



**CARLO ERBA**  
*Reagents operates with  
 a Certified Quality  
 Management System*

