



## Tapones para tubos de ensayo desechables

**Tapones para tubos de ensayo desechables**  
Poliétileno.



Code	Description	Packaging
LLG09231641	Caps 15-17 mm for test tubesneutral/100	100 pz.
LLG09231642	Caps 11-13 mm for test tubesneutral/100	100 pz.
LLG09231643	Caps 15-17 mm for test tubesblue, P/100	100 pz.
LLG09231644	Caps 11-13 mm for test tubesblue, P/100	100 pz.
LLG09231645	Caps 15-17 mm for test tubesyellow,/100	100 pz.
LLG09231646	Caps 11-13 mm for test tubesyellow,/100	100 pz.
LLG09231647	Caps 15-17 mm for test tubesred, PE/100	100 pz.
LLG09231648	Caps 11-13 mm for test tubesred, PE/100	100 pz.
LLG09231649	Caps 15-17 mm for test tubesgreen, /100	100 pz.
LLG09231650	Caps 11-13 mm for test tubesgreen, /100	100 pz.



**CARLO ERBA Reagents SA**  
 Calle Filadors 35,  
 6º Planta 5º Puerta  
 08208 Sabadell (BCN)  
 Tel. +34 93 693 37 35  
[www.carloerbareagents.com](http://www.carloerbareagents.com)



**CARLO ERBA**  
*Reagents operates with  
 a Certified Quality  
 Management System*

