



Tapones para tubos de ensayo desechables

Tapones para tubos de ensayo desechables
Poliétileno.



Code	Description	Packaging
LLG09231641	Caps 15-17 mm for test tubesneutral/100	100 pz.
LLG09231642	Caps 11-13 mm for test tubesneutral/100	100 pz.
LLG09231643	Caps 15-17 mm for test tubesblue, P/100	100 pz.
LLG09231644	Caps 11-13 mm for test tubesblue, P/100	100 pz.
LLG09231645	Caps 15-17 mm for test tubesyellow,/100	100 pz.
LLG09231646	Caps 11-13 mm for test tubesyellow,/100	100 pz.
LLG09231647	Caps 15-17 mm for test tubesred, PE/100	100 pz.
LLG09231648	Caps 11-13 mm for test tubesred, PE/100	100 pz.
LLG09231649	Caps 15-17 mm for test tubesgreen, /100	100 pz.
LLG09231650	Caps 11-13 mm for test tubesgreen, /100	100 pz.



CARLO ERBA Reagents SA
 Calle Filadors 35,
 6º Planta 5º Puerta
 08208 Sabadell (BCN)
 Tel. +34 93 693 37 35
www.carloerbareagents.com



CARLO ERBA
*Reagents operates with
 a Certified Quality
 Management System*

