

**Termohigrómetro****TFA Dostmann****Termohigrómetro**

Pelo sintético, carcasa con anillo de bronce.

Rango de medición

Humedad relativa del aire: 0 hasta 100 %

Temperatura: -15 hasta 55 °C

Pesa: 65 g

**CARLO ERBA**
REAGENTS

Code	Description	Packaging
LLG09726403	Thermohygrometers, analogue0-100%, -5..	1 pz



CARLO ERBA Reagents SA
Calle Filadors 35,
6º Planta 5º Puerta
08208 Sabadell (BCN)
Tel. +34 93 693 37 35
www.carloerbareagents.com



CARLO ERBA
*Reagents operates with
a Certified Quality
Management System*

