

**Termohigrómetro "Klimatherm"****TFA Dostmann****Termohigrómetro "Klimatherm"**

Carcasa con anillo de bronce.

**Rango de medición****Humedad relativa del aire: 15 hasta 100 %****Temperatura: 0 hasta 40 °C****Pesa: 61 g****CARLO ERBA**  
REAGENTS

Code	Description	Packaging
LLG09726404	Thermohygrometers Klimatherm15-100%, 0.	1 pz



**CARLO ERBA Reagents SA**  
Calle Filadors 35,  
6º Planta 5º Puerta  
08208 Sabadell (BCN)  
Tel. +34 93 693 37 35  
[www.carloerbareagents.com](http://www.carloerbareagents.com)



**CARLO ERBA**  
*Reagents operates with  
a Certified Quality  
Management System*

