

**Termohigrómetro "Klimatherm"****TFA Dostmann****Termohigrómetro "Klimatherm"**

Carcasa con anillo de bronce.

Rango de medición

Humedad relativa del aire: 15 hasta 100 %

Temperatura: 0 hasta 40 °C

Pesa: 61 g



Code	Description	Packaging
LLG09726404	Thermohygrometers Klimatherm15-100%, 0.	1 pz



CARLO ERBA Reagents SA
Calle Filadors 35,
6º Planta 5º Puerta
08208 Sabadell (BCN)
Tel. +34 93 693 37 35
www.carloerbareagents.com



CARLO ERBA
*Reagents operates with
a Certified Quality
Management System*

