



**Tubos de ensayo, vidrio sodocálcico**

**Tubos de ensayo, vidrio sodocálcico**



Code	Description	Packaging
LLG04658073	Test Tubes 25x7,00x0,6-0,7mm Soda /1000	1000 pz.
LLG04658074	Test Tubes 35x6,50x0,6-0,7mm Soda /1000	1000 pz.
LLG04658075	Test Tubes 35x8,00x0,6-0,7mm Soda /1000	1000 pz.
LLG04658076	Test Tubes 40x8.00x0.6-0.7mm Soda l/200	200 pz.
LLG04658077	Test Tubes 40x9.50x0.6-0.7mm Soda l/500	500 pz.
LLG04658078	Test Tubes 50x6,50x0,6-0,7mm Soda /1000	1000 pz.
LLG04658079	Test Tubes 50x9,50x0,7-0,8mm Soda l/200	200 pz.
LLG04658080	Test Tubes 55x11,00x0,7-0,8mm Soda /500	500 pz.
LLG04658081	Test Tubes 70x11,00x0,7-0,8mm Soda /200	200 pz.
LLG04658082	Test Tubes 75x12,00x0,5-0,6mm Soda /250	250 pz.
LLG04658083	Test Tubes 100x12,00x0,5-0,6mm Soda/250	250 pz.
LLG04658084	Test Tubes 100x13,00x0,5-0,6mm Soda/250	250 pz.
LLG04658085	Test Tubes 120x12,00x0,6-0,7mm Soda/500	500 pz.
LLG04658086	Test Tubes 125x16,00x0,6-0,7mm Soda/250	250 pz.
LLG04658088	Test Tubes 150x18,00x0,8-0,9mm Soda/250	250 pz.
LLG04658089	Test Tubes 160x16,00x0,6-0,7mm Soda/250	250 pz.
LLG04658090	Test Tubes 180x18,00x0,7-0,8mm Soda/100	100 pz.
LLG04658091	Test Tubes 120x16,00x0,5-0,6mm Soda/100	100 pz.
LLG04658092	Test Tubes 160x16,00x0,5-0,6mm Soda/100	100 pz.
LLG04658093	Test Tubes 180x18,00x0,6-0,7mm Soda/100	100 pz.



Code	Description	Packaging
LLG04658094	Test Tubes 180x20,00x0,7-0,8mm Soda/100	100 pz.
LLG06232384	Test tubes 75 x 11.50 x 0.7-0.8 mm /250	250 pz.
LLG06267191	Test tubes,soda glass,without rim,1/250	250 pz.



**CARLO ERBA Reagents SA**  
 Calle Filadors 35,  
 6º Planta 5º Puerta  
 08208 Sabadell (BCN)  
 Tel. +34 93 693 37 35  
[www.carloerbareagents.com](http://www.carloerbareagents.com)



**CARLO ERBA**  
*Reagents operates with  
 a Certified Quality  
 Management System*

