



## Tubos de ensayo, vidrio de borosilicato 5.1

Tubos de ensayo, vidrio de borosilicato 5.1



| Code        | Description                             | Packaging |
|-------------|-----------------------------------------|-----------|
| LLG04658095 | Test Tubes 75x10,00x0,5-0,6mm Boro /250 | 250 pz.   |
| LLG04658096 | Test Tubes 125x16,00x0,6-0,7mm Boro/250 | 250 pz.   |
| LLG04658097 | Test Tubes 150x16,00x0,6-0,7mm Boro/250 | 250 pz.   |
| LLG04658099 | Test Tubes 70x8,00x0,4-0,5mm Boro 5/100 | 100 pz.   |
| LLG04658100 | Test Tubes 75x10,00x0,4-0,5mm Boro /100 | 100 pz.   |
| LLG04658101 | Test Tubes 75x12,00x0,4-0,5mm Boro /100 | 100 pz.   |
| LLG04658102 | Test Tubes 100x10,00x0,4-0,5mm Boro/100 | 100 pz.   |
| LLG04658103 | Test Tubes 100x12,00x0,4-0,5mm Boro/100 | 100 pz.   |
| LLG04658104 | Test Tubes 100x16,00x0,5-0,6mm Boro/100 | 100 pz.   |
| LLG04658105 | Test Tubes 130x14,00x0,4-0,5mm Boro/100 | 100 pz.   |
| LLG04658106 | Test Tubes 130x16,00x0,5-0,6mm Boro/100 | 100 pz.   |
| LLG04658107 | Test Tubes 150x20,00x0,5-0,6mm Boro/100 | 100 pz.   |
| LLG04658108 | Test Tubes 150x25,00x0,6-0,7mm Boro /50 | 50 pz.    |
| LLG04658109 | Test Tubes 180x18,00x0,5-0,6mm Boro/100 | 100 pz.   |
| LLG04658110 | Test Tubes 180x20,00x0,5-0,6mm Boro/100 | 100 pz.   |
| LLG04658111 | Test Tubes 200x25,00x0,6-0,7mm Boro /50 | 50 pz.    |
| LLG04658112 | Test Tubes 200x30,00x0,7-0,8mm Boro /50 | 50 pz.    |
| LLG06237590 | Test tubes 150x20mm, wall 0.8-0.9 m/100 | 100 pz.   |
| LLG06265921 | Test tubes 100x16mm, wall 0,6-0,7 m/250 | 250 pz.   |
| LLG07623956 | Reagent tubes 12 x 75mm, wall 0.5-0/250 | 250 pz.   |





**CARLO ERBA Reagents SA**  
Calle Filadors 35,  
6º Planta 5º Puerta  
08208 Sabadell (BCN)  
Tel. +34 93 693 37 35  
[www.carloerbareagents.com](http://www.carloerbareagents.com)



**CARLO ERBA**  
*Reagents operates with  
a Certified Quality  
Management System*

