



**Tapones de rosca para frascos de boca ancha, PP/LDPE**

**Tapones de rosca para frascos de boca ancha, PP/LDPE**  
Tapón de PP, junta de LDPE.



Code	Description	Packaging
LLG06272535	Screw-cap DIN 45, black PPN, LDPE insert	1 pz
LLG09072164	Srew cap PP for 30ml wide mouth bottles	1 pz.
LLG09072165	Screw cap PP for 50,100 and 125 mlwide	1 pz.
LLG09072166	Screw cap PP for 250 and 500ml bottlesD	1 pz.
LLG09072167	Screw cap PP for bottles 1000ml/DIN 68	1 pz.



**CARLO ERBA Reagents SA**  
 Calle Filadors 35,  
 6º Planta 5º Puerta  
 08208 Sabadell (BCN)  
 Tel. +34 93 693 37 35  
[www.carloerbareagents.com](http://www.carloerbareagents.com)



**CARLO ERBA**  
*Reagents operates with  
 a Certified Quality  
 Management System*

