



Tapones de rosca para frascos de boca ancha, PP/LDPE

Tapones de rosca para frascos de boca ancha, PP/LDPE

Tapón de PP, junta de LDPE.



Code	Description	Packaging
LLG06272535	Screw-cap DIN 45, black PPN, LDPE insert	1 pz
LLG09072164	Srew cap PP for 30ml wide mouth bottles	1 pz.
LLG09072165	Screw cap PP for 50,100 and 125 mlwide	1 pz.
LLG09072166	Screw cap PP for 250 and 500ml bottlesD	1 pz.
LLG09072167	Screw cap PP for bottles 1000ml/DIN 68	1 pz.



CARLO ERBA Reagents SA
 Calle Filadors 35,
 6º Planta 5º Puerta
 08208 Sabadell (BCN)
 Tel. +34 93 693 37 35
www.carloerbareagents.com



CARLO ERBA
*Reagents operates with
 a Certified Quality
 Management System*

