



Bidones industriales, HDPE

Kautex

Bidones industriales, HDPE

Bidones de 2 a 30 litros de capacidad.

Pedir los cierres por separado.



Code	Description	Packaging
LLG09140021	Industrial jerrycan without closure2 l,	1 pz.
LLG09140022	Industrial jerrycan without closure2,5	1 pz.
LLG09140023	Industrial jerrycan without closure3 l,	1 pz.
LLG09140025	Industrial jerrycan without closure5 l,	1 pz.
LLG09140026	Industrial jerrycan without closure6 l,	1 pz.
LLG09140027	Industrial jerrycan without closure8, l	1 pz.
LLG09140030	Industrial jerrycan without closure10 l	1 pz.
LLG09140031	Industrial jerrycan without closure25 l	1 pz.
LLG09140032	Industrial jerrycan without closure20 l	1 pz.
LLG09140033	Industrial jerrycan without closure30 l	1 pz.
LLG09140037	Industrial jerrycans without closure20	1 pz.



CARLO ERBA Reagents SA
 Calle Filadors 35,
 6º Planta 5º Puerta
 08208 Sabadell (BCN)
 Tel. +34 93 693 37 35
www.carloerbareagents.com



CARLO ERBA
*Reagents operates with
 a Certified Quality
 Management System*

