



**Accesorios para tubos de PEEK**

Acoplamiento hermético para tubos de HPLC fabricados en PEEK.



Code	Description	Packaging
LLG04664157	PEEK Fingertight fittings 1/16 to 1	10 pz.
LLG04664158	PEEK Fingertight fittings 1/32 to 1	10 pz.
LLG04664160	PEEK Fingertight fittings 1/32 to 6	5 pz.



**CARLO ERBA Reagents SA**  
 Calle Filadors 35,  
 6ª Planta 5ª Puerta  
 08208 Sabadell (BCN)  
 Tel. +34 93 693 37 35  
[www.carloerbareagents.com](http://www.carloerbareagents.com)



**CARLO ERBA**  
*Reagents operates with  
 a Certified Quality  
 Management System*

