



Accesorios para tubos de PEEK

Acoplamiento hermético para tubos de HPLC fabricados en PEEK.



Code	Description	Packaging
LLG04664157	PEEK Fingertight fittings 1/16 to 1	10 pz.
LLG04664158	PEEK Fingertight fittings 1/32 to 1	10 pz.
LLG04664160	PEEK Fingertight fittings 1/32 to 6	5 pz.



CARLO ERBA Reagents SA
 Calle Filadors 35,
 6º Planta 5º Puerta
 08208 Sabadell (BCN)
 Tel. +34 93 693 37 35
www.carloerbareagents.com



CARLO ERBA
*Reagents operates with
 a Certified Quality
 Management System*

