



Accesorio para tintómetro óptico, Lovibond<sup>®</sup> modelo F

Tintometer

Accesorio para tintómetro óptico, Lovibond<sup>®</sup> modelo F



Code	Description	Packaging
LLG09947346	Replacement sample chamberinsert for /3	3 pz.
LLG09947347	Accessory pack for solidsample measurem	1 pz.



**CARLO ERBA Reagents SA**  
Calle Filadors 35,  
6º Planta 5º Puerta  
08208 Sabadell (BCN)  
Tel. +34 93 693 37 35  
[www.carloerbareagents.com](http://www.carloerbareagents.com)



**CARLO ERBA**  
*Reagents operates with  
a Certified Quality  
Management System*

