

**Vaso y bolas de molienda para mini-molino PULVERISETTE 23****Fritsch****Vaso y bolas de molienda para mini-molino PULVERISETTE 23****Vaso de molienda**|5 ml|10 ml|15 ml**Capacidad útil (volumen de la muestra)**|0,1-1ml|0,2-1 ml|0,5-5 ml**Diámetro de las bolas**|||

5 mm|20|30|60

10 mm|1|3|5

15 mm|-|1|2



Code	Description	Packaging
LLG09739343	Grinding bowls 5 ml volume PTFE	1 pz.
LLG09739348	Grinding bowl 15ml hardened steel/1	1 pz.
LLG09739349	Grinding bowl 10ml hardened steel/1	1 pz.



**CARLO ERBA Reagents SA**  
 Calle Filadors 35,  
 6º Planta 5º Puerta  
 08208 Sabadell (BCN)  
 Tel. +34 93 693 37 35  
[www.carloerbareagents.com](http://www.carloerbareagents.com)



**CARLO ERBA**  
*Reagents operates with  
 a Certified Quality  
 Management System*

