

**Vaso y bolas de molienda para mini-molino PULVERISETTE 23****Fritsch****Vaso y bolas de molienda para mini-molino PULVERISETTE 23****Vaso de molienda** | 5 ml | 10 ml | 15 ml**Capacidad útil (volumen de la muestra)** | 0,1-1 ml | 0,2-1 ml | 0,5-5 ml**Diámetro de las bolas** | | |

5 mm | 20 | 30 | 60

10 mm | 1 | 3 | 5

15 mm | - | 1 | 2



Code	Description	Packaging
LLG09739343	Grinding bowls 5 ml volume PTFE	1 pz.
LLG09739348	Grinding bowl 15ml hardened steel/1	1 pz.
LLG09739349	Grinding bowl 10ml hardened steel/1	1 pz.



CARLO ERBA Reagents SA
 Calle Filadors 35,
 6º Planta 5º Puerta
 08208 Sabadell (BCN)
 Tel. +34 93 693 37 35
www.carloerbareagents.com



CARLO ERBA
*Reagents operates with
 a Certified Quality
 Management System*

