



Sonotrodo para procesador ultrasonido compacto UP200S

Hielscher Ultrasonics

Sonotrodo para procesador ultrasonido compacto UP200S



Code	Description	Packaging
LLG09106607	Sonotrode S1	1 pz.
LLG09106608	Sonotrode S3	1 pz.
LLG09106609	Sonotrode S7	1 pz.
LLG09106610	Sonotrode S14	1 pz.
LLG09106611	Sonotrode S14Dwith glass flow vessel GD	1 pz.



CARLO ERBA Reagents SA
Calle Filadors 35,
6º Planta 5º Puerta
08208 Sabadell (BCN)
Tel. +34 93 693 37 35
www.carloerbareagents.com



CARLO ERBA
*Reagents operates with
a Certified Quality
Management System*

