



Sonotrodo para procesador ultrasonido compacto UP200S

Hielscher Ultrasonics

Sonotrodo para procesador ultrasonido compacto UP200S



| Code | Description | Packaging |
|-------------|---|-----------|
| LLG09106607 | Sonotrode S1 | 1 pz. |
| LLG09106608 | Sonotrode S3 | 1 pz. |
| LLG09106609 | Sonotrode S7 | 1 pz. |
| LLG09106610 | Sonotrode S14 | 1 pz. |
| LLG09106611 | Sonotrode S14Dwith glass flow vessel GD | 1 pz. |



CARLO ERBA Reagents SA
 Calle Filadors 35,
 6º Planta 5º Puerta
 08208 Sabadell (BCN)
 Tel. +34 93 693 37 35
www.carloerbareagents.com



CARLO ERBA
*Reagents operates with
 a Certified Quality
 Management System*

