



Silicone rubber stoppers

Silicone rubber stoppers

Transparent. To DIN 12871. Solid. Operating temperature up to 220 °C.
Supplied singly.



Code	Description	Packaging
LLG09230503	Silikonstopfen 3,5x6,5x15mmacc.to DIN 1	1 pz
LLG09230505	Silicone rubber stoppers, 5 x 9 x 20 mm	1 pz
LLG09230508	Silicone rubber stoppers, 8 x 12 x 20 m	1 pz
LLG09230510	Silicone rubber stoppers, 10,5 x 14,5 x	1 pz
LLG09230512	Silicone rubber stoppers, 12,5 x 16,5 x	1 pz
LLG09230514	Silicone rubber stoppersm 14 x 18 x 20	1 pz
LLG09230517	Silicone stopper 17x22x25mmacc. to DIN	1 pz
LLG09230518	Silicone rubber stoppers, 18 x 24 x 30	1 pz
LLG09230521	Silicone rubber stoppers, 21 x 27 x 30	1 pz
LLG09230523	Silicone rubber stoppers, 23 x 29 x 30	1 pz
LLG09230526	Silicone rubber stoppers, 26 x 32 x 30	1 pz
LLG09230529	Silicone rubber stoppers, 29 x 35 x 30	1 pz
LLG09230531	Silicone rubber stoppers, 31 x 38 x 35	1 pz
LLG09230536	Silicone rubber stoppers, 36 x 44 x 40	1 pz
LLG09230541	Silicone rubber stoppers, 41 x 49 x 40	1 pz
LLG09230547	Silicone rubber stopper, 47 x 55 x 40 m	1 pz
LLG09230550	Silicone rubber stoppers, Bottom diam. 50.5 mm, Top diam. 59.5 mm, Height 45 mm	1 pz
LLG09230556	Silicone rubber stoppers, Bottom diam. 56 mm, Top diam. 65 mm, Height 45 mm	1 pz
LLG09230560	Silicone rubber stoppers, Bottom diam. 60 mm, Top diam. 70 mm, Height 50 mm	1 pz
LLG09230564	Silicone rubber stoppers, Bottom diam. 64.5 mm, Top diam. 75.5 mm, Height 55 mm	1 pz



Code	Description	Packaging
LLG09230571	Silicone rubber stoppers, Bottom diam. 71 mm, Top diam. 83 mm, Height 60 mm	1 pz
LLG09230579	Silicone rubber stoppers, Bottom diam. 79 mm, Top diam. 92 mm, Height 65 mm	1 pz



CARLO ERBA Reagents S.A.S.
 Chaussée du Vexin – 27106 Val
 de Reuil cedex
 N° TVA: FR 63391048824
 Tél : +33 (0)2 32 09 20 00
www.carloerbareagents.com



CARLO ERBA
*Reagents operates with
 a Certified Quality
 Management System*

