

**Tapones de silicona**

Deutsch &amp; Neumann

**Tapones de silicona**

Transparentes. Según DIN 12871. Sólidos. Temperatura de trabajo hasta 220 °C. Se suministran individualmente.



Code	Description	Packaging
LLG09230503	Silikonstopfen 3,5x6,5x15mmacc.to DIN 1	1 pz.
LLG09230505	Silicone rubber stoppers, 5 x 9 x 20 mm	1 pz.
LLG09230508	Silicone rubber stoppers, 8 x 12 x 20 m	1 pz.
LLG09230510	Silicone rubber stoppers, 10,5 x 14,5 x	1 pz.
LLG09230512	Silicone rubber stoppers, 12,5 x 16,5 x	1 pz.
LLG09230514	Silicone rubber stoppersm 14 x 18 x 20	1 pz.
LLG09230517	Silicone stopper 17x22x25mmacc. to DIN	1 pz.
LLG09230518	Silicone rubber stoppers, 18 x 24 x 30	1 pz.
LLG09230521	Silicone rubber stoppers, 21 x 27 x 30	1 pz.
LLG09230523	Silicone rubber stoppers, 23 x 29 x 30	1 pz.
LLG09230526	Silicone rubber stoppers, 26 x 32 x 30	1 pz.
LLG09230529	Silicone rubber stoppers, 29 x 35 x 30	1 pz.
LLG09230531	Silicone rubber stoppers, 31 x 38 x 35	1 pz.
LLG09230536	Silicone rubber stoppers, 36 x 44 x 40	1 pz.
LLG09230541	Silicone rubber stoppers, 41 x 49 x 40	1 pz.
LLG09230547	Silicone rubber stopper, 47 x 55 x 40 m	1 pz.
LLG09230550	Silicone rubber stoppers, 50,5 x 59,5 x	1 pz.
LLG09230556	Silicone rubber stoppers, 56 x 65 x 45	1 pz.
LLG09230560	Silicone rubber stoppers, 60 x 70 x 50	1 pz.
LLG09230564	Silicone rubber stoppers, 64,5 x 75,5 x	1 pz.



Code	Description	Packaging
LLG09230571	Silicone rubber stoppers, 71 x 83 x 60	1 pz.
LLG09230579	Silicone rubber stoppers79x92x65mm	1 pz.



**CARLO ERBA Reagents SA**  
Calle Filadors 35,  
6º Planta 5º Puerta  
08208 Sabadell (BCN)  
Tel. +34 93 693 37 35  
[www.carloerbareagents.com](http://www.carloerbareagents.com)



**CARLO ERBA**  
*Reagents operates with  
a Certified Quality  
Management System*

