



Recipientes para molino de bolas de alta energía Emax

Recipientes para molino de bolas de alta energía Emax



Code	Description	Packaging
LLG09738776	Grinding beaker 50ml, stainless steel f	1 pz.
LLG09738777	Grinding beaker 125ml, stainless steel	1 pz.
LLG09738778	Gridning beaker 50ml, zirconium oxide f	1 pz.
LLG09738779	Grinding beaker 125ml, zirconium oxide	1 pz.
LLG09738780	Grinding beaker 50ml, tungsten carbide	1 pz.



CARLO ERBA Reagents SA
 Calle Filadors 35,
 6º Planta 5ª Puerta
 08208 Sabadell (BCN)
 Tel. +34 93 693 37 35
www.carloerbareagents.com



CARLO ERBA
*Reagents operates with
 a Certified Quality
 Management System*

