



Recipientes para molino de bolas de alta energía Emax

Recipientes para molino de bolas de alta energía Emax



| Code | Description | Packaging |
|-------------|---|-----------|
| LLG09738776 | Grinding beaker 50ml, stainless steel f | 1 pz. |
| LLG09738777 | Grinding beaker 125ml, stainless steel | 1 pz. |
| LLG09738778 | Gridning beaker 50ml, zirconium oxide f | 1 pz. |
| LLG09738779 | Grinding beaker 125ml, zirconium oxide | 1 pz. |
| LLG09738780 | Grinding beaker 50ml, tungsten carbide | 1 pz. |



CARLO ERBA Reagents SA
 Calle Filadors 35,
 6º Planta 5º Puerta
 08208 Sabadell (BCN)
 Tel. +34 93 693 37 35
www.carloerbareagents.com



CARLO ERBA
*Reagents operates with
 a Certified Quality
 Management System*

