

**Cápsulas de cristalización, DURAN<sup>®</sup>****DURAN Group****Cápsulas de cristalización, DURAN<sup>®</sup>**

- DIN 12338
- Con pico
- Cápsulas a partir de 100 ml con espacio para rotulación
- Tipo de vidrio I/vidrio neutro según USP, EP y JP
- Esterilizables en autoclave



Code	Description	Packaging
LLG09144024	Crystallizing dishes,DURAN,with spout,c	1 pz
LLG09144032	Crystallizing dishes,DURAN,with spout,c	1 pz
LLG09144034	Crystallizing dishes,DURAN,with spout,c	1 pz
LLG09144038	Crystallizing dishes,DURAN,with spout,c	1 pz
LLG09144041	Crystallizing dishes,DURAN,with spout,c	1 pz
LLG09144044	Crystallizing dishes,DURAN,with spout,c	1 pz
LLG09144049	Crystallizing dishes,DURAN,with spout,c	1 pz
LLG09144054	Crystallizing dishes,DURAN,with spout,c	1 pz
LLG09144059	Crystallizing dishes,DURAN,with spout,c	1 pz
LLG09144063	Crystallizing dishes,DURAN,with spout,c	1 pz



**CARLO ERBA Reagents SA**  
 Calle Filadors 35,  
 6º Planta 5º Puerta  
 08208 Sabadell (BCN)  
 Tel. +34 93 693 37 35  
[www.carloerbareagents.com](http://www.carloerbareagents.com)



**CARLO ERBA**  
*Reagents operates with  
 a Certified Quality  
 Management System*

