

**Cápsulas de cristalización, DURAN<sup>®</sup>****DURAN Group****Cápsulas de cristalización, DURAN<sup>®</sup>**

- DIN 12338
- Con pico
- Cápsulas a partir de 100 ml con espacio para rotulación
- Tipo de vidrio I/vidrio neutro según USP, EP y JP
- Esterilizables en autoclave



| Code        | Description                             | Packaging |
|-------------|-----------------------------------------|-----------|
| LLG09144024 | Crystallizing dishes,DURAN,with spout,c | 1 pz      |
| LLG09144032 | Crystallizing dishes,DURAN,with spout,c | 1 pz      |
| LLG09144034 | Crystallizing dishes,DURAN,with spout,c | 1 pz      |
| LLG09144038 | Crystallizing dishes,DURAN,with spout,c | 1 pz      |
| LLG09144041 | Crystallizing dishes,DURAN,with spout,c | 1 pz      |
| LLG09144044 | Crystallizing dishes,DURAN,with spout,c | 1 pz      |
| LLG09144049 | Crystallizing dishes,DURAN,with spout,c | 1 pz      |
| LLG09144054 | Crystallizing dishes,DURAN,with spout,c | 1 pz      |
| LLG09144059 | Crystallizing dishes,DURAN,with spout,c | 1 pz      |
| LLG09144063 | Crystallizing dishes,DURAN,with spout,c | 1 pz      |



**CARLO ERBA Reagents SA**  
Calle Filadors 35,  
6º Planta 5ª Puerta  
08208 Sabadell (BCN)  
Tel. +34 93 693 37 35  
[www.carloerbareagents.com](http://www.carloerbareagents.com)



**CARLO ERBA**  
*Reagents operates with  
a Certified Quality  
Management System*

