

**Cápsulas de cristalización, DURAN[®]****DURAN Group****Cápsulas de cristalización, DURAN[®]**

- DIN 12338
- Con pico
- Cápsulas a partir de 100 ml con espacio para rotulación
- Tipo de vidrio I/vidrio neutro según USP, EP y JP
- Esterilizables en autoclave



Code	Description	Packaging
LLG09144024	Crystallizing dishes,DURAN,with spout,c	1 pz
LLG09144032	Crystallizing dishes,DURAN,with spout,c	1 pz
LLG09144034	Crystallizing dishes,DURAN,with spout,c	1 pz
LLG09144038	Crystallizing dishes,DURAN,with spout,c	1 pz
LLG09144041	Crystallizing dishes,DURAN,with spout,c	1 pz
LLG09144044	Crystallizing dishes,DURAN,with spout,c	1 pz
LLG09144049	Crystallizing dishes,DURAN,with spout,c	1 pz
LLG09144054	Crystallizing dishes,DURAN,with spout,c	1 pz
LLG09144059	Crystallizing dishes,DURAN,with spout,c	1 pz
LLG09144063	Crystallizing dishes,DURAN,with spout,c	1 pz



CARLO ERBA Reagents SA
Calle Filadors 35,
6º Planta 5º Puerta
08208 Sabadell (BCN)
Tel. +34 93 693 37 35
www.carloerbareagents.com



CARLO ERBA
*Reagents operates with
a Certified Quality
Management System*

