



Accesorios para autoclaves, serie HG

HMC-EUROPE

Accesorios para autoclaves, serie HG



Code	Description	Packaging
LLG06251590	Printer HG series	1 pz.
LLG09842420	Flexible temperature sensor	1 pz.
LLG09842421	automatic water feedfor HG 50	1 pz.
LLG09842422	Exhaust air filter for HG seriefor the	1 pz.
LLG09842423	Filling carrier HG 50/80Ø 345 x 225 mm	1 pz.
LLG09842424	Filling carrierwith drip trayØ 345 x 22	1 pz.
LLG09842425	Filling bucketØ 345 x 195 mm H	1 pz.



CARLO ERBA Reagents SA
 Calle Filadors 35,
 6º Planta 5º Puerta
 08208 Sabadell (BCN)
 Tel. +34 93 693 37 35
www.carloerbareagents.com



CARLO ERBA
*Reagents operates with
 a Certified Quality
 Management System*

