



## Viscosímetros, sin boquilla

### Viscosímetros, sin boquilla

Para boquillas intercambiables de 2, 3, 4, 5, 6, 7 y 8 mm de salida. Pedir por separado. Viscosímetros y boquillas similares a DIN 53211.



| Code        | Description                             | Packaging |
|-------------|-----------------------------------------|-----------|
| LLG09149860 | Viscosity cup,DIN 53211,Aluminium,witho | 1 pz.     |
| LLG09149872 | Orifices for viscosity cups,DIN 53211,  | 1 pz.     |
| LLG09149873 | Orifices for viscosity cups,DIN 53211,  | 1 pz.     |
| LLG09149874 | Orifices for viscosity cups,DIN 53211,  | 1 pz.     |
| LLG09149875 | Orifices for viscosity cups,DIN 53211,  | 1 pz.     |
| LLG09149876 | Orifices for viscosity cups,DIN 53211,  | 1 pz.     |
| LLG09149877 | Orifices for viscosity cups,DIN 53211,  | 1 pz.     |
| LLG09149878 | Orifices for viscosity cups,DIN 53211,  | 1 pz.     |



**CARLO ERBA Reagents SA**  
Calle Filadors 35,  
6º Planta 5º Puerta  
08208 Sabadell (BCN)  
Tel. +34 93 693 37 35  
[www.carloerbareagents.com](http://www.carloerbareagents.com)



**CARLO ERBA**  
*Reagents operates with  
a Certified Quality  
Management System*

