



Viscosímetros, sin boquilla

Viscosímetros, sin boquilla

Para boquillas intercambiables de 2, 3, 4, 5, 6, 7 y 8 mm de salida. Pedir por separado.
Viscosímetros y boquillas similares a DIN 53211.



Code	Description	Packaging
LLG09149860	Viscosity cup,DIN 53211,Aluminium,witho	1 pz.
LLG09149872	Orifices for viscosity cups,DIN 53211,	1 pz.
LLG09149873	Orifices for viscosity cups,DIN 53211,	1 pz.
LLG09149874	Orifices for viscosity cups,DIN 53211,	1 pz.
LLG09149875	Orifices for viscosity cups,DIN 53211,	1 pz.
LLG09149876	Orifices for viscosity cups,DIN 53211,	1 pz.
LLG09149877	Orifices for viscosity cups,DIN 53211,	1 pz.
LLG09149878	Orifices for viscosity cups,DIN 53211,	1 pz.



CARLO ERBA Reagents SA
 Calle Filadors 35,
 6º Planta 5º Puerta
 08208 Sabadell (BCN)
 Tel. +34 93 693 37 35
www.carloerbareagents.com



CARLO ERBA
*Reagents operates with
 a Certified Quality
 Management System*

