




Soluciones estándar para conductímetros

Valores indicados a 25°C, Precisión: +/- 2%.



Code	Description	Packaging
LLG09041053	Conductivity test solution LF 9903x6 am	18 pz.
LLG09041055	Conductivity test solution LF 9923x6 am	18 pz.
LLG09041056	Conductivity test solution LF 1024 K250	1 pz.



CARLO ERBA Reagents SA
Calle Filadors 35,
6º Planta 5º Puerta
08208 Sabadell (BCN)
Tel. +34 93 693 37 35
www.carloerbareagents.com



CARLO ERBA
*Reagents operates with
a Certified Quality
Management System*

