



## Bandejas para instrumental, resina de melamina

### Bandejas para instrumental, resina de melamina

Blancas, de resina de melamina. Autoclavables (121°C).



| Code        | Description                             | Packaging |
|-------------|-----------------------------------------|-----------|
| LLG07055383 | Instrument tray 30x19x2cmmelamine resin | 1 pz.     |
| LLG07608352 | Instrument tray 19x15x2cmmelamine resin | 1 pz.     |
| LLG09200324 | Instrument tray, 24x18x2 cmmelamine res | 1 pz.     |
| LLG09200327 | Instrument tray, 27x21x2 cmmelamine res | 1 pz.     |
| LLG09200328 | Instrument tray 30x15x2cmmelamine resin | 1 pz.     |
| LLG09200335 | Instrument tray, 36x24x2 cmmelamine res | 1 pz.     |
| LLG09200345 | Instrument tray, 43x29x2 cmmelamine res | 1 pz.     |



**CARLO ERBA Reagents SA**  
Calle Filadors 35,  
6º Planta 5ª Puerta  
08208 Sabadell (BCN)  
Tel. +34 93 693 37 35  
[www.carloerbareagents.com](http://www.carloerbareagents.com)



**CARLO ERBA**  
*Reagents operates with  
a Certified Quality  
Management System*

