



Probetas graduadas, PP, forma alta, clase B, graduación moldeada

VITLAB

Probetas graduadas, PP, forma alta, clase B, graduación moldeada

Forma alta. Muy transparentes. Con una escala en relieve y marcas de anillo en los puntos primarios de la escala, calibradas en "In". Tolerancias de clase B según DIN 12681 / ISO 6706. La base hexagonal con tacos inferiores proporciona una alta estabilidad. La tensión térmica hasta 80°C no provoca que se superen permanentemente los límites de tolerancia. Producto apto para uso alimentario según el reglamento (UE) N° 10/2011.



Code	Description	Packaging
LLG07075031	Measuring cylinder 10 mlPP, tall form,c	1 pz.
LLG09274420	Measuring cylinder 25 mlPP, tall form,c	1 pz.
LLG09274428	Measuring cylinder 50 mlPP, tall form,c	1 pz.
LLG09274438	Measuring cylinder 100 mlPP, tall form,	1 pz.
LLG09274448	Measuring cylinder 250 mlPP, tall form,	1 pz.
LLG09274454	Measuring cylinder 500 mlPP, tall form,	1 pz.
LLG09274462	Measuring cylinder 1000 mlPP, tall form	1 pz.
LLG09274464	Measuring cylinder 2000 mlPP, tall form	1 pz.



CARLO ERBA Reagents SA
 Calle Filadors 35,
 6º Planta 5º Puerta
 08208 Sabadell (BCN)
 Tel. +34 93 693 37 35
www.carloerbareagents.com



CARLO ERBA
*Reagents operates with
 a Certified Quality
 Management System*

