




## Cucharas para fósforo

### Cucharas para fósforo

Cápsula de acero inoxidable. 16 mm de Ø. Longitud del mango, 450 mm.



Code	Description	Packaging
LLG09150610	Deflagrating spoon,steel,diam. 16 mm,/1	1 pz.



**CARLO ERBA Reagents SA**  
Calle Filadors 35,  
6º Planta 5º Puerta  
08208 Sabadell (BCN)  
Tel. +34 93 693 37 35  
[www.carloerbareagents.com](http://www.carloerbareagents.com)



**CARLO ERBA**  
*Reagents operates with  
a Certified Quality  
Management System*

