



Accesorios para bomba para barril GLP 25

Grün-Pumpen

Accesorios para bomba para barril GLP 25



| Code        | Description                              | Packaging |
|-------------|--|-----------|
| LLG09001400 | Motor for barrel pump GLP25230 V         | 1 pz      |
| LLG09001401 | Pump station barrel pump GLP25length 50  | 1 pz      |
| LLG09001402 | Pump station barrel pump GLP25length 70  | 1 pz      |
| LLG09001403 | Pump station barrel pump GLP25length 10  | 1 pz      |
| LLG09001404 | "angled tube PPAG 3/4\"", length 250 mm" | 1 pz      |
| LLG09001405 | Stopper PPD= 65-70 mm                    | 1 pz      |
| LLG09001406 | Pump nozzle, PPwith 1.5 m PVC tube, DN 1 | 1 pz      |



**CARLO ERBA Reagents SA**  
 Calle Filadors 35,  
 6ª Planta 5ª Puerta  
 08208 Sabadell (BCN)  
 Tel. +34 93 693 37 35  
[www.carloerbareagents.com](http://www.carloerbareagents.com)



**CARLO ERBA**  
*Reagents operates with  
 a Certified Quality  
 Management System*

