



**Cuentahilos de precisión metálico**

**Eschenbach**

**Cuentahilos de precisión metálico**

- Marco de latón, cromado, acabado mate seda
- Lentes de vidrio de sílice, biconvexas



| Code        | Description                             | Packaging |
|-------------|-----------------------------------------|-----------|
| LLG09151423 | Thread counters,brass,height 23 mm,fiel | 1 pz      |
| LLG09151427 | Thread counters,brass,height 27 mm,fiel | 1 pz      |
| LLG09151433 | Thread counters,brass,height 33 mm,fiel | 1 pz      |
| LLG09151442 | Thread counters,brass,height 42 mm,fiel | 1 pz      |
| LLG09151453 | Thread counters,brass,height 53 mm,fiel | 1 pz      |



**CARLO ERBA Reagents SA**  
Calle Filadors 35,  
6º Planta 5ª Puerta  
08208 Sabadell (BCN)  
Tel. +34 93 693 37 35  
[www.carloerbareagents.com](http://www.carloerbareagents.com)



**CARLO ERBA**  
*Reagents operates with  
a Certified Quality  
Management System*

