



Cuentahilos de precisión metálico

Eschenbach

Cuentahilos de precisión metálico

- Marco de latón, cromado, acabado mate seda
- Lentes de vidrio de sílice, biconvexas



Code	Description	Packaging
LLG09151423	Thread counters,brass,height 23 mm,fiel	1 pz
LLG09151427	Thread counters,brass,height 27 mm,fiel	1 pz
LLG09151433	Thread counters,brass,height 33 mm,fiel	1 pz
LLG09151442	Thread counters,brass,height 42 mm,fiel	1 pz
LLG09151453	Thread counters,brass,height 53 mm,fiel	1 pz



CARLO ERBA Reagents SA
 Calle Filadors 35,
 6º Planta 5º Puerta
 08208 Sabadell (BCN)
 Tel. +34 93 693 37 35
www.carloerbareagents.com



CARLO ERBA
*Reagents operates with
 a Certified Quality
 Management System*

