



## Filtro de aspiración para disolventes

### Filtro de aspiración para disolventes

Filtro de disolventes de HPLC para tubos de 3,2 mm de diámetro externo.



Code	Description	Packaging
LLG04005890	"HPLC solvent filters for 1/8" (3.2mm)/	5 pz.
LLG04005891	HPLC solvent filters for 1/8 (3.2mm	5 pz.



**CARLO ERBA Reagents SA**  
 Calle Filadors 35,  
 6º Planta 5º Puerta  
 08208 Sabadell (BCN)  
 Tel. +34 93 693 37 35  
[www.carloerbareagents.com](http://www.carloerbareagents.com)



**CARLO ERBA**  
*Reagents operates with  
 a Certified Quality  
 Management System*

