



Filtro de aspiración para disolventes

Filtro de aspiración para disolventes

Filtro de disolventes de HPLC para tubos de 3,2 mm de diámetro externo.



Code	Description	Packaging
LLG04005890	"HPLC solvent filters for 1/8" (3.2mm)/	5 pz.
LLG04005891	HPLC solvent filters for 1/8 (3.2mm	5 pz.



CARLO ERBA Reagents SA
Calle Filadors 35,
6º Planta 5º Puerta
08208 Sabadell (BCN)
Tel. +34 93 693 37 35
www.carloerbareagents.com



CARLO ERBA
*Reagents operates with
a Certified Quality
Management System*

