



Probetas de mezcla, vidrio de borosilicato 3.3, forma alta, clase A, graduación azul

Isolab

Probetas de mezcla, vidrio de borosilicato 3.3, forma alta, clase A, graduación azul

- Certificado de conformidad
- DIN EN ISO 4788
- Con tapón de PE
- Pie hexagonal
- Las inscripciones están impresas con esmalte azul
- El certificado de lote se puede descargar desde el sitio web



| Code | Description | Packaging |
|-------------|---|-----------|
| LLG09275070 | Mixing cylinder 10 ml, tall formglass, | 1 pz. |
| LLG09275071 | Mixing cylinder 25 ml, tall formglass, | 1 pz. |
| LLG09275072 | Mixing cylinder 50 ml, tall formglass, | 1 pz. |
| LLG09275073 | Mixing cylinder 100 ml, tall formglass, | 1 pz. |
| LLG09275074 | Mixing cylinder 250 ml, tall formglass, | 1 pz. |
| LLG09275075 | Mixing cylinder 500 ml, tall formglass, | 1 pz. |
| LLG09275076 | Mixing cylinder 1000 ml, tall formglass | 1 pz. |



CARLO ERBA Reagents SA
 Calle Filadors 35,
 6º Planta 5ª Puerta
 08208 Sabadell (BCN)
 Tel. +34 93 693 37 35
www.carloerbareagents.com



CARLO ERBA
*Reagents operates with
 a Certified Quality
 Management System*

