



Probetas de mezcla, vidrio de borosilicato 3.3, forma alta, clase A, graduación azul

Isolab

Probetas de mezcla, vidrio de borosilicato 3.3, forma alta, clase A, graduación azul

- Certificado de conformidad
- DIN EN ISO 4788
- Con tapón de PE
- Pie hexagonal
- Las inscripciones están impresas con esmalte azul
- El certificado de lote se puede descargar desde el sitio web



Code	Description	Packaging
LLG09275070	Mixing cylinder 10 ml, tall formglass,	1 pz.
LLG09275071	Mixing cylinder 25 ml, tall formglass,	1 pz.
LLG09275072	Mixing cylinder 50 ml, tall formglass,	1 pz.
LLG09275073	Mixing cylinder 100 ml, tall formglass,	1 pz.
LLG09275074	Mixing cylinder 250 ml, tall formglass,	1 pz.
LLG09275075	Mixing cylinder 500 ml, tall formglass,	1 pz.
LLG09275076	Mixing cylinder 1000 ml, tall formglass	1 pz.



CARLO ERBA Reagents SA
 Calle Filadors 35,
 6º Planta 5ª Puerta
 08208 Sabadell (BCN)
 Tel. +34 93 693 37 35
www.carloerbareagents.com



CARLO ERBA
*Reagents operates with
 a Certified Quality
 Management System*

