



Accesorios para TITRONIC[®] 500 y TitroLine[®] 7000 / 7500 KF / 7750 / 7800

SI Analytics

Accesorios para TITRONIC[®] 500 y TitroLine[®] 7000 / 7500 KF / 7750 / 7800



Code	Description	Packaging
LLG09863154	Magnetic stirrer TM 235115-230 V	1 pz.
LLG09863155	Change unit 5ml WA 05complete	1 pz.
LLG09863156	Change unit 10ml WA 10complete	1 pz.
LLG09863157	Change unit 20ml WA 20complete	1 pz.
LLG09863158	Change unit 50ml WA 50complete	1 pz.
LLG09863159	USB keyboard TZ 3835	1 pz.



CARLO ERBA Reagents SA
 Calle Filadors 35,
 6º Planta 5º Puerta
 08208 Sabadell (BCN)
 Tel. +34 93 693 37 35
www.carloerbareagents.com



CARLO ERBA
*Reagents operates with
 a Certified Quality
 Management System*

